


IMPLANT EXTRACTION SYSTEM



bti.[®]

Biotechnology
Institute
Human Technology



Human Technology

IMPLANT EXTRACTION SYSTEM

KIT KEXIM



THE ONLY SYSTEM on the market that guarantees atraumatic extraction of internal and external connection implants, in any situation.

THE FIRST SYSTEM ON THE MARKET PATENTED BY BTI FOR SAFE EXTRACTION

Extractor kit

Five exclusive extractors indicated for the extraction of osseointegrated implants when they are no longer functional or are unnecessary in a new prosthetic situation.

Trephine drill kit

An atraumatic trephine drill kit for an ultrafine cut to release the area with maximum osseointegration, preserving the surrounding bone, only drilling the first 3 millimetres.

EXTRACTOR KIT PARTS

COUNTER-TORQUE WRENCH,
200 NCM.

Wrench with large torque that breaks the join between the bone and impact, making the extraction process easier. When torque over 200 Ncm is applied, the wrench forms an angle and prevents progress. Back to the original position using a cylindrical component in the kit.



MDPT

Handle to carry the extractor once the extender has been fitted.



EXTENDERS

Set of counter-torque wrench extenders (long, short, extralong) for different clinical situations.



EXTRACTORS*

Five parts with varying widths and lengths to adapt to the range of implants on the market.

If necessary due to the depth of the connection, they can be cut to size using a diamond disc cutter.


They are disposable to ensure maximum performance.

* Never exceed 70 Ncm with the narrowest one (EXIME - N° 0). In this case, replace LLT200 with LLMQ.

AN EASY AND PRECISE PROCEDURE


1

1. Choose the appropriate implant extractor. Choose the right extractor for the shape of the internal and external implant connections.



2

2. Manually insert the extractor in the implant thread, using the handle-adapter, ensuring it fits in an anti-clockwise direction.



3

3. Connect the counter-torque wrench and turn it anti-clockwise.



4

4. Apply torque gently in an anti-clockwise direction until the implant is fully extracted.



If the torque of the counter-torque wrench (200 Ncm.) is not sufficient for the extraction, Re-set it and re-attempt a maximum of three times, maintaining the torque for a few seconds before releasing. If the implant can still not be extracted, drill around it using a trephine drill to a depth of just 2-3mm and repeat the wrench implant extraction process.

TREPHINE DRILL KIT

BTI SOLUTION FOR ATRAUMATIC EXTRACTION

TREPHINE DRILLS

Cylindrical and hollow drills to cut the bone surrounding the implant. They are delivered in four different diameters, depending on the width of the implants:

- **Tiny:** FTEX37-5
- **Universal:** FTEX43-5
- **Universal plus:** FTEX51-5
- **Wide:** FTEX57-5




ADVANTAGES:

· They allow an ultrafine cut along the first 5 millimetres because they are narrower at that level, preserving the surrounding bone. This is particularly important in the case of resorbed maxilla.

TREPHINE DRILLING TO THE MILLIMETRE

Trephine drilling to just 2-3mm reduces maximum torque to less than half, to break the implant-bone connection with guaranteed extraction success.





BTI Biotechnology Institute
San Antonio, 15 · 5º
01005 Vitoria-Gasteiz
(Álava) · SPAIN
Tel: +34 945 140 024
Fax: +34 945 135 203
pedidos@bticomercial.com

www.bti-biotechnologyinstitute.com

USA
1730 Walton Road
Suite 110
Blue Bell, PA 19422-1802 · USA
Tel: (12) 156 464 067
Fax: (12) 156 464 066
info@bti-implant.us

UK
870 The Crescent
Colchester Business Park · Colchester
Essex CO49YQ · United Kingdom
Tel: (44) 01206580160
Fax: (44) 01206580161
info@bti-implant.co.uk

GERMANY
Mannheimer Str. 17
75179 Pforzheim · Germany
Tel: (49) 7231 428 06-0
Fax: (49) 7231 428 06-15
info@bti-implant.de

ITALY
Piazzale Piola, 1
20131 Milan · Italy
Tel: (39) 02 7060 5067
Fax: (39) 02 7063 9876
bti.italia@bti-implant.it

MEXICO
Ejercito Nacional Mexicano 351, 3A
Col. Granada Delegación Miguel Hidalgo
México DF · CP 11520 · México
Tel: (52) 55 52502964
Fax: (52) 55 55319327
bti.mexico@bti-implant.com

PORTUGAL
Praça Mouzinho de Albuquerque 113, 5º
4100-359 Porto · Portugal
Tel: (351) 22 120 1373
Fax: (351) 22 120 1311
bti.portugal@bticomercial.com



PRODUCT NEWS
BIOTECHNOLOGY INSTITUTE

SCAN THIS CODE with your mobile phone to access all the information of **BTI Biotechnology Institute**.

